

+PATIENT DATA

+NO ONE GETS YOU CLOSER

## TrialMax Slate™ Site-Based ePRO

### Overview

At CRF Health we understand that the data collection needs of clinical studies are varied and complex. The demand for electronic methods for the collection of Patient Reported Outcomes (PRO) data has increased, and so too have the questionnaire requirements. As a result, larger sized screens are required for certain site-based PRO instruments. We have extended our TrialMax® and TrialStudio® to support large touch-screen Tablet PCs. TrialMax Slate's larger screens allow more complex questionnaires to be completed with ease and with minimal changes to the original instrument layout. The flexibility of stylus and touch screen allow for more natural information entry, simple navigation as well as brilliant color and clarity for easier reading and usability.



### What is TrialMax Slate™?

TrialMax Slate™ is intended to be used at the study sites during subject visits for completion of QoL and other PRO instruments during those visits. TrialMax Slate™ is part of the TrialMax® product family and can be deployed independently or in combination with TrialMax Touch™ home-based eDiaries.

### Benefits of TrialMax Slate™:

- **Usability:** the large display area allows for multiple questions on one screen, making it easier to complete and avoids having to change the original instrument layout
- **Design:** The TrialMax® eDiary design tool, TrialStudio® offers a truly collaborative approach to the eDiary design process with real-time development capabilities
- **Data Management:** The data captured by TrialMax Slate™ is always clean because of real-time edit checks and the source data is electronic, which eliminates the need for source data verification
- **Real-time Data Access:** Sites, CRAs and the study team can access all captured data in real time, improving the information flow within the study
- **Large Screen Support:** TrialMax Slate™ allows the ePRO instrument to be designed with large buttons. This increases usability for patients who may be visually impaired or have muscle or joint problems
- **Graphical Interface:** Allows for the incorporation of text, graphics and image maps - making the questionnaire more appealing and easy-to-use



+PATIENT DATA

+NO ONE GETS YOU CLOSER

## TrialMax Slate™ Site-Based ePRO

### TrialMax Slate™ Technical Details

#### Data Sending Options with TrialMax Slate™ ePRO Devices

- Wi-Fi
- Cellular data / 3G wireless connections
- Landline modem
- LAN port

#### TrialMax Slate™ Displays

- Ranges from 10.1" to 12.1"
- Multiple resolutions and aspect ratios available
- Finger or stylus operated touch screens

#### Battery Life

Battery life varies by model, usage characteristics and parameter settings, but all devices are commercial laptop computers and can be used unplugged for several hours before requiring recharging. In most cases, the devices will be plugged into an outlet at the site, but can easily be used independently to provide patient privacy during assessment completion.

#### Language Support

CRF Health's TrialMax® ePRO technology has been used in more than 70 countries, on six continents and 68 regional languages, including several Indian dialects. The TrialMax® system offers virtually unlimited language support.

#### Memory and Performance



CRF Health thoroughly evaluates all devices to ensure that supported models perform well in clinical studies. Our TrialMax Slate™ devices are responsive and easy-to-use, with a robust design that holds up to use in clinical trials.

Modern Tablet PCs have ample memory and fast processors, ensuring the devices can handle complex calculations and the needs of your protocol.

#### Suitable PRO Instruments

The TrialMax Slate™ system is suitable for all site-based PROs as well as ClinROs (Clinician Reported Outcomes). CRF Health has deployed a large number of different instruments for a variety of indications ranging from Alzheimer's disease to Oncology. TrialMax Slate™ supports a large number of distinct users and roles, including patients, caregivers and clinicians as required by today's more complex protocols.

CRF Health can assist you in determining whether TrialMax Slate™ is a good option for your clinical trial.

#### Integration with EDC and Medical Devices

- Medidata Rave®
- TrialMax Synapse™
- Open Interface
- Glucometer
- Spirometer

